

	Technical Data Sheet
Use in	<ul style="list-style-type: none"> Pharmaceutical Industry Medical Device Industry Cosmetic Industry For industrial, laboratory & research applications only Medium according to EP and USP
Use for	<ul style="list-style-type: none"> Washing solution
Typical composition per liter	<p>NaCl 9 g</p> <p>This medium can be adjusted / or supplemented according to the performance criteria required.</p>
Filling volume	<ul style="list-style-type: none"> 1000 mL
Bottle format	<ul style="list-style-type: none"> 1000 mL laboratory glass screw cap bottle Type I glass (borosilicate glass) Bottle opening about 30 mm Colour of cap: rot GL45 screw cap with 3 integrated septa
Bottles per tray	<ul style="list-style-type: none"> 6 bottles on a plastic tray wrapped with shrink foil
Shelf life	<ul style="list-style-type: none"> 18 months from production date
Storage conditions	<ul style="list-style-type: none"> Recommended storage temperature: 2 - 25 °C Should be stored at temperatures as stable as possible Store protected from light exposure
Label	<ul style="list-style-type: none"> On the side Contain autoclave indicator
Label information	<ul style="list-style-type: none"> Product name: Physiol.NaCl 1000 mL Expiry date: YYYYMMDD → MMM in letters (e.g.: 2023Nov04) Lot-number Individual number Barcode
Barcode	<ul style="list-style-type: none"> 2-dimensional (data matrix), 20 digits: Digits 1-3: Art.-No. Digits 4-9: Lot-Number Digits 10-14: Individual-Number Digits 15-20: Date (YYMMDD)
Delivery	<ul style="list-style-type: none"> Temperature controlled delivery on request For shipments of larger amounts plastic pallets in Euro-size can be used

Technical Data Sheet	
Bottle information	<ul style="list-style-type: none"> Label contains autoclaving indicator (brown → green) Bottles are incubated at 25 – 35 °C for at least 48 hours after autoclaving and then packed Bottles are not touched any more by hand after autoclaving
Place of production	PharmaMedia Dr. Müller GmbH Gustav-Throm-Str. 1, 69181 Leimen - Germany

Quality control, Certificates									
Certificates	<p>Each lot of product can be obtained with a certificate of analysis (CoA):</p> <table border="1"> <tr> <th colspan="2">Physico-chemical test parameters:</th></tr> <tr> <td>Filling volume</td><td>980 – 1040 mL</td></tr> <tr> <td colspan="2">Sterility control</td></tr> <tr> <td colspan="2">≥ 7 days at 53-55 °C with bioindicators, no growth</td></tr> </table>	Physico-chemical test parameters:		Filling volume	980 – 1040 mL	Sterility control		≥ 7 days at 53-55 °C with bioindicators, no growth	
Physico-chemical test parameters:									
Filling volume	980 – 1040 mL								
Sterility control									
≥ 7 days at 53-55 °C with bioindicators, no growth									
Release of negative pressure in media bottles	During the autoclaving process, chemical reactions inside the bottles may result in a slight vacuum. Please assure that the vacuum is released without contaminating the bottle. Ideally, the vacuum is released by puncturing the septum with an aeration needle equipped with a sterile filter prior to opening a bottle.								
Certificate of origin	<p>All media lots produced by PMM can be obtained with a Certificate of Origin (CoO). All animal derived raw materials are specified as follows:</p> <ul style="list-style-type: none"> Raw material Tissue Animal source Country of origin Infectivity category (acc. to TSE guideline: EMA/410/01 current version) 								
BSE policy	<ul style="list-style-type: none"> In compliance with the current note for guidance on minimizing the risk of transmitting animal spongiform encephalopathy via human or veterinary medicinal products, we check the CoO of raw material in respect to the specified animal source, the country of origin and the infectivity category. We neither store or process ruminant raw materials obtained from high infectivity tissues (IA) nor ruminant raw materials whose animal source originates from countries or regions with an undetermined risk (cat C/GBR IV). 								

	Safety Data
Toxic ingredients	<ul style="list-style-type: none">• None
Basic composition	<ul style="list-style-type: none">• See typical composition
Solvent content	<ul style="list-style-type: none">• None
Safety data sheet required	<ul style="list-style-type: none">• Not mandatorily required